

Vampire-Dracula Facelift Study (iFRAX Micro- Fractional RF System + Mitate C3 Mask)

**Objective: Administer series of
three iFRAX procedures in
combination with Mitate C3 Mask
and provide photographs and
clinical feedback.**

Study Objective and Treatment Schedule

- Treatment on subject whose skin is dry, dehydrated, dull, and “loose”. Also presenting freckles and scars on right side of face
- 1st iFRAX Treatment Date: May 30, 2013
- Apply Mitate C3 mask: June 10, 2013
- 2nd iFRAX Treatment Date: June 19~23, 2013
- Apply Mitate C3 mask: July 1 & July 10, 2013
- 3rd iFRAX Treatment Date: July 18~20, 2013
- Apply Mitate C3 mask: July 30, 2013

Application of Topical Anesthetic Cream

EMLA 5% lidocaine prilocaine



Topical Is Wiped Off & Ready For Treatment



Treatment Parameter Settings for “White Zone”

- Interval (Duration) - 5%
- Energy Level - 20%
- Needle Depth – 1.5



iFRAX Treatment of “White Zone”



iFRAX Treatment of “White Zone” Observations

- The subject reported feeling no pain during treatment
- Forehead area appeared more red than other areas
- Only very minor bleeding apparent on forehead area
- Redness on forehead reduced quickly after treatment was finished

After Treating “White Zone” with iFRAX



iFRAX Treatment Settings for “Pink” Zone

- Interval (Duration) – 5%
- Energy Level – 25%
- Needle Depth – 3 mm



Treating “Pink Zone” with iFRAX System



After Treatment of “Pink Zone” with iFRAX



- Note: Performed a second pass over at uneven scar area on the right side face (as shown)

Application of Mitate C3 Mask



Sonophoresis with SkinKarma Ultratoner



After Mitate C3 Mask and Sonophoresis Treatment

- During Mitate C3 mask treatment the subject felt tingling on her face (this feeling lasted for several minutes)
- During SkinKarma Ultratoner treatment, the subject reported no unusual feeling or discomfort
- After treatment, the subject felt no discomfort or “tingling”
- The subject reported that her face felt cool and tighter after the removal of the mask

Immediately after Treatment



Immediately after Treatment

- The outer corner of the subject's eyes and forehead appeared full, firm and hydrated
- The subject reported feeling a firming and tightening effect on her face
- The subject's skin appeared brighter and tighter

24 Hours after Treatment



24 Hours after Treatment

- OB liposome and sun protection applied at home after treatment
- The subject's skin felt and appeared more hydrated and calm

48 Hours after Treatment



48 Hours after Treatment

- The subject continued to use OB liposome and hydrating mask at home after the treatment
- The subject's skin appeared tighter, brighter, and more hydrated after the treatment
- The subject's pore size diminished and she was pleased with the results of her treatment